

THE K⁺, Na⁺ TRANSPORT ATP-ase

The semipermeable plasma membrane of eukaryotic cells serves as a permeability barrier separating the cells from their *surroundings*. The flow of molecules and ions between the cell and its environment is precisely regulated by specific transport systems. In most animal cells the concentration of K⁺ is high, and the concentration of Na⁺ relative to the external medium is low. These ionic gradients are generated and maintained against passive diffusion by the Na/K pump (K⁺, Na⁺ transport ATP-ase). The energy needed to keep the pump going is provided by the hydrolysis of ATP. The transport can be inhibited by ouabain - a cardiotonic steroid - that binds to the external surface of the membrane and inhibits the dephosphorylation of the K⁺, Na⁺ transport ATP-ase. The ATP-ase activity is determined by assaying the amount of liberated inorganic phosphate (P_i) using the Fiske-Subbarow colour reaction.

Solutions:

- 1.) 0.5 M TRIS HCl buffer pH=7.4
- 2.) 0.05 M MgCl₂
- 3.) mixture of 1 M NaCl and of 0.2 M KCl solution
- 4.) 1 M NaCl
- 5.) 0.2 M KCl
- 6.) 5 mM Ouabain (Strophantine)
- 7.) 5 mM ATP
- 8.) 20 % TCA
- 9.) 2.5 % Ammonium molybdate
- 10.) 1 % ascorbic acid
- 11.) microsomal fraction of grey matter of porcine brain as enzyme source

Preparation of microsomal fraction of porcine brain

Isolating solution: 30 mM TRIS-HCl
0,25 M Saccharose pH= 7.4

Homogenate (20 %) of grey matter of porcine brain was prepared in **an** isolating solution. Homogenate was centrifuged at 6,000. RPM for 30 min. at 4°C (Janetzky K26 centrifuge) and the pellet was discarded. Supernatant was recentrifuged at 40,000 RPM for 60 min at 4°C (Beckman Ultracentrifuge). The pellet consisting of the microsomal fraction was resuspended in the isolating solution and its protein content was determined with the Biuret method. Microsomal fraction was diluted to a 2 mg/ml of protein concentration in isolating solution.

Experiment

Pipette into test tubes the following ingredients:

Materials/#	1	2	3	4	5	6	7
TRIS-HCl ml	0.25	0.25	0.25	0.25	0.25	0.25	0.25
MgCl ₂ ml	0.1	0.1	0.1	0.1	0.1	0.1	0.1
NaCl-KCl ml	-	-	0.1	-	0.1	-	-
NaCl ml	-	-	-	-	-	0.1	-
KCl ml	-	-	-	-	-	-	0.1
Oubain ml	-	-	-	0.2	0.2	-	-
H ₂ O ml	0.3	0.3	0.2	0.1	-	0.2	0.2
ATP ml	0.25	0.25	0.25	0.25	0.25	0.25	0.25
Microsomal suspensions ml	-	0.1	0.1	0.1	0.1	0.1	0.1

Incubate the tubes for 20 minutes at 37 °C .

Add the followings to the test tubes:

TCA ml	0.5	0.5	0.5	0.5	0.5	0.5	0.5
Micros. susp.	0.1	-	-	-	-	-	-

Transfer the content of the tubes into Eppendorf tubes and centrifuge with 10 000 rpm for 2 min.

Determination of the amount of liberated phosphate (P_i):

Add into 7 numbered test tubes the followings:

Supernatants (from the corresponding Eppendorf tubes)	0.3 ml
H ₂ O	3,2 ml
Ammonium molybdate	1,0 ml
Ascorbic acid	0.5 ml

After incubation at 37 °C for 5 min, cool to room temperature.

Measure the absorption of samples at 720 nm. Water serves as blank. Using the calibration curve calculate the amounts of P_i.

Calculate the percentage of the activity derived from the K⁺/Na⁺ transport ATP-ase compared to the total ATP-ase activity.